

Monitoring of HTS Compound Library Quality via a High-Resolution Image Acquisition and Processing Instrument

Pierre Baillargeon,¹ Louis Scampavia,^{1,2} Ross Einsteder,¹ and Peter Hodder^{1,2*}

¹*Scripps Research Institute Molecular Screening Center, Lead Identification, Translational Research Institute, Jupiter, FL*

²*Department of Molecular Therapeutics, Scripps Florida, Jupiter, FL*

Keywords:
compound management, HTS library, machine vision, precipitate detection, colored compound, volume detection

This report presents the high-resolution image acquisition and processing instrument for compound management applications (HI-API-CM). The HI-API-CM combines imaging spectroscopy and machine-vision analysis to perform rapid assessment of *high-throughput screening* (HTS) compound library quality. It has been customized to detect and classify typical artifacts found in HTS compound library microtiter plates (MTPs). These artifacts include (1) insufficient volume of liquid compound sample, (2) compound precipitation, and (3) colored compounds that interfere with HTS assay detection format readout. The HI-API-CM is also configured to automatically query and compare its analysis results to data stored in a LIMS or corporate database, aiding in the detection of compound registration errors. To demonstrate its capabilities, several compound plates ($n = 5760$ wells total) containing different artifacts were measured via automated HI-API-CM analysis, and results compared with those obtained by manual (visual) inspection. In all cases, the instrument demonstrated high

fidelity (99.8% empty wells; 100.1% filled wells; 94.4% for partially filled wells; 94.0% for wells containing colored compounds), and in the case of precipitate detection, the HI-API-CM results significantly exceeded the fidelity of visual observations (220.0%). As described, the HI-API-CM allows for noninvasive, nondestructive MTP assessment with a diagnostic throughput of about 1 min per plate, reducing analytical expenses and improving the quality and stewardship of HTS compound libraries. (JALA 2011;16:197–203)

INTRODUCTION

With the evolution of reliable *high-throughput screening* (HTS) automation,¹ large libraries of drug-like compounds are screened for pharmacologic activity in pharmaceutical, biotechnology, and academic screening facilities throughout the world. Because the cost of a typical compound averages between \$10 and \$100/mg,² many HTS centers select microtiter plates (MTPs) as the preferred vehicle for compound storage and handling. Although compound library size varies³ at different research facilities, the vast majority are solvated in dimethyl sulfoxide (DMSO). As DMSO-solvated liquids, compounds are easily formatted into high-density MTPs (e.g., 384-well, 1536-well), which are ideally suited for containing the minuscule amount (10–100 μ L) of compound necessary for

*Correspondence: Peter Hodder, Ph.D., Lead Identification Division, Translational Research Institute, Scripps Florida, Jupiter, FL 33458; Phone: +1.561.228.2100; Fax: +1.561.228.3054; E-mail: hodderp@scripps.edu

2211-0682/\$36.00

Copyright © 2011 by the Society for Laboratory Automation and Screening

doi:10.1016/j.jala.2011.02.004

HTS studies. In many cases, a research facility may have several copies of its library, either stored directly on the HTS platform or in longer-term cold storage.^{3,4}

As HTS libraries are solvated in DMSO, they are subject to the introduction of artifacts (Fig. 1). Unfortunately, DMSO is hygroscopic and susceptible to water assimilation. Therefore, compounds stored as DMSO-solvated HTS libraries degrade over time through exposure, such as freeze-thaw cycles or storage in humid environments.^{5–8} As a consequence, artifact may be introduced to an HTS library, in the form of sample precipitation or compound hydrolysis. Additional artifacts may be introduced by typical activities in a compound management operation. For example, a liquid handler may fail to correctly aliquot a compound from a source MTP to a destination plate; this error may result in an empty or partially filled well in the destination plate. Alternatively, a particular compound sample may be over-sampled by a hit-picker, resulting in an empty well in the source MTP. A colored compound may also be considered as an artifact; for example, its absorbance spectrum may overlap with the fluorescence or luminescence generated by an HTS assay reagent. Lastly, archival errors in an outdated or poorly maintained laboratory information management system (LIMS) can also complicate HTS screening efforts from the generation of mismatched biological and chemical data.

The challenges in quality assurance (QA) and quality control (QC) of compound libraries have led many facilities to use traditional analytical instruments for periodic “spot”

inspections of an HTS library during its lifecycle. Unfortunately, these technologies are ill suited for screening the vast number of compounds contained in HTS libraries. At present, liquid chromatography–mass spectrometry (LC–MS) instrumentation is the gold standard that combines the physical separation capabilities of liquid chromatography (i.e., purity) with the mass analysis capabilities of MS (i.e., compound identification). In practice, LC–MS is limited to small batch validations, because it requires extensive analysis times, destructive sampling, and large consumable expenses.⁹ Nephelometry is limited to precipitate detection and can be unreliable.^{10,11} Acoustic auditing devices are used with some success to monitor solvate composition and volume in MTPs.^{12–15} However, this type of instrumentation cannot distinguish different well artifacts (i.e., precipitate/crystals from oil, bubbles, or color) and is limited to analysis of proprietary MTPs, requiring costly replating of large HTS libraries to realize its application. Moreover, the sample throughput of all these technologies is too slow to be capable of performing rapid, routine surveys of entire library collections. To date, the most popular method to diagnose well artifacts is through visual inspection. Obviously, this manual technique is fatiguing and subjective. In larger HTS libraries, visual inspection is virtually impossible.⁵

In an effort to address this unmet need, we have developed a novel machine-vision-based technology referred to as the high-resolution image acquisition and processing instrument for compound management applications (HI-API-CM). The HI-API-CM provides an innovative solution to MTP




Figure 1. Common artifacts found in MTP-based HTS compound libraries. Artifacts observed in individual wells of a 384-well MTP are illustrated: (A) an empty well, (B) a well containing colored compound, (C) a well with bubbles, (D) a partially filled well, (E) and a well containing precipitate. For comparison, an image of a “normal” DMSO-filled well (F) is included. The HI-API-CM detects and classifies all these artifacts in an automated fashion with high fidelity.

inspection through the use of dual-image analysis. By incorporating visible (VIS) and infrared (IR) imaging spectroscopy, it can be used for rapid detection of common artifacts found in MTP-based HTS compound libraries. The instrument performs three simultaneous measurements in a noninvasive fashion, namely: (1) discrimination of empty and partially filled MTP wells, (2) detection of insoluble compound crystals/precipitates, and (3) identification of potential colorimetric interferents. It can display results through a user interface (UI), or compare its results to those in a LIMS for fully automated QA methods. In this article, a description of the HI-API-CM, its software and performance in a pilot study is presented.

INSTRUMENT OVERVIEW

HI-API-CM uses dual-image analyses, using VIS and IR images taken with a charge-coupled device (CCD) camera. Fabricated and integrated at our research facility, it consists of a color CCD camera (detector), automated MTP loader, telecentric imaging lens, plate barcode reader, and dual-mode illuminator. Instrument control and image analysis computers are located underneath the optical/mechanical tower within a mobile rack enclosure (Fig. 2). This unit is compact and portable, readily carted to serve multiple sites as needed.

A distinguishing feature of the HI-API-CM is its use of a custom telecentric lens coupled to a high-resolution,

progressive-scanning five-megapixel CCD camera. This optical arrangement is capable of imaging the entire MTP without a need for repositioning. With a depth of field of ~ 12 mm, the lens produces negligible perspective and focal errors, and the instrument is able to resolve features to ~ 60 μm with minimal vignetting. Consequently, the HI-API-CM is capable of detecting artifacts missed by manual inspections.

A customized light source was designed, fabricated, and integrated on the HI-API-CM to provide illumination using either VIS or IR light. As seen in Figure 3, IR imaging reveals hidden artifacts (i.e., precipitate/crystals) that may otherwise be masked by color. The illuminator requires no warm-up period for even illumination and instantaneously toggles between VIS and IR modes. This feature allows for multiple image analyses across the MTP field in ~ 1 min, and is not limited to any appreciable extent by MTP format density (i.e., 96-, 384-, 1536-wells). Transmission mode illumination is used in the presented setup, with images taken from the bottom of the MTP.

Of equal importance to the HI-API-CM hardware is its use of customized *machine-vision* software for image processing, and a UI for operating the instrument and reviewing the results. When applied to the analysis of HTS compound libraries, machine-vision methods facilitate the detection of features that evade more conventional detection methods. To illustrate this utility, Figure 4 presents VIS light snapshots of different MTP wells from an HTS compound




Figure 2. The HI-API-CM. Shown on the left is the HI-API-CM with a fully functioning plate loader, barcode reader, and servers in the mobile computer rack below. On the right, a CAD rendering of the prototype's interior with a light illumination configured for a "bottom plate read" mode is illustrated. A telecentric lens allows for a full plate imaging by the detector. Complete analysis is approximately a minute per plate.




Figure 3. HI-API-CM distinguishes colored compounds from precipitate. All images from 384-well MTPs. VIS-light images (A, C) of colored compounds can mask the presence of precipitate. However, the HI-API-CM's IR images of these same wells clearly identifies the top well (images A and B) as containing insoluble precipitate, whereas the bottom well (images C and D) as containing solubilized compound.

library, and intermediate machine-vision processing steps that result in the automated classification of the sample contained in each well. Machine-vision processing begins as the MTP image is digitally partitioned into individual wells, each with a predefined “area of interest” (AOI) to be analyzed. After selection of the AOI, the first step in analysis is the extraction of color information, which is spectroscopically resolved in terms of pixel intensity at each color channel. After the extraction of color information from the AOI, image subtraction is performed from an archived control plate, consisting of empty wells. This background correction facilitates further analysis of the image based only on the well contents found in the AOI. Grayscale profiles can be created that transect pixel variances across the well, which can be further simplified to binarized images. The software can convert these binarized images into digital “object signatures.” Predefined digital thresholds are then applied by the software to distinguish the different well types. For example, an empty compound well will result in a black binarized image with no object signature. If the well is not empty, it will yield a different result at each stage of image processing, and be automatically classified as appropriate.

UI SOFTWARE

A UI and associated computer hardware facilitates the seamless deployment of the HI-API-CM. Two computers are located on the HI-API-CM platform (Fig. 2). The first is a “client” that controls the instrument and runs the image acquisition software, whereas the other hosts the Apache web and structured query language database (MySQL)

applications, which are used for the storage and retrieval of archived images and metadata. As unprocessed raw images are acquired within $\sim 250\ \mu\text{m}$, most of the image analysis time ($\sim 1\ \text{min}$) is expended for machine-vision processes and database storage.

Figure 5 presents actual MTP images from the customized HI-API-CM software, as viewed through the UI. A 384-well polypropylene plate (Greiner, USA) containing DMSO-solvated compounds with various artifacts is shown in the unprocessed raw image. This raw image is juxtaposed to images where the HI-API-CM has displayed the results of its analysis. As seen in the figure, classification of wells containing colored samples, precipitate, and different volumes of sample are easily viewed through the interactive UI software and easily identified by the user. The HI-API-CM is also capable of comparing its analysis results to those archived in a LIMS database via an SQL query. Similar to the flagging of other artifacts, the UI software highlights wells that have conflicting status between the HI-API-CM and LIMS databases. Alternatively, the user can elect an ad hoc analysis where no database comparison is performed. If the LIMS comparison is not chosen, any well that deviates from those properties of a “normal” well (i.e., a well filled with DMSO) will be classified and tagged by the HI-API-CM.

Postanalysis, the UI software permits manual corrections, which facilitates software development and training efforts. However, the HI-API-CM has been designed to run fully automated with a capability of exporting flat-text files that summarize results. The HI-API-CM database archive also allows users to search and select images from any particular plate (retrievable by barcode or HI-API-CM assay date) as might be necessary for long-term QC/QA inspection and intervention.

HI-API-CM PILOT STUDY

To demonstrate “proof-of-concept,” a small inventory of 15 HTS compound library plates (384-well format) with known artifacts were analyzed using the HI-API-CM. The HI-API-CM results were compared with manual inspections (performed visually by eye). The HI-API-CM took 15 min to complete its analysis, including plate loading and unloading; in contrast, the manual inspections took 75 min. The HI-API-CM successfully detected many of the artifacts found through manual inspection with high fidelity (Table 1).

Discussion of the HI-API-CM's detection fidelity provides insight into different aspects of the HI-API-CM software and hardware design. For example, the HI-API-CM had excellent fidelity in classifying empty and full wells; automation of partially filled well classification exhibited slightly lower fidelity. The disagreement between manual and automated analysis can be attributed to the HI-API-CM platform's software design. Specifically, when the HI-API-CM software is unable to confidently classify an artifact, it defaults to a conservative assignment of well status. Although infrequent, the effects of this classification scheme are most apparent when the volume




Figure 4. Machine-vision techniques resolve common artifacts found in HTS compound libraries. Examples of a “normal” (DMSO-filled) well and four different artifact types found in individual wells of a 384-well MTP are presented. From left to right, the wells are a DMSO-filled well, a well partially filled with DMSO, a well containing colored compound, a well with precipitated compound, and an empty well. All wells can be analyzed and categorized via the HI-API-CM’s machine-vision analysis technology.

of sample in an interrogated well is on the boundary between the HI-API-CM’s predefined “empty” and “partially filled” assignment cutoffs. In the case of the MTPs used for this pilot study, a threshold volume of at least 2 μL was found to be necessary for 100% accuracy in the HI-API-CM’s partially filled well determination. Before executing the pilot study, the HI-API-CM software was programmed to assign an “empty” status to any well whose volume was between 0 and 2 μL . This programming is reflected in the 94.4% agreement between manual and automated classification and to a lesser extent in the fidelity of other volumetric classification results. As the development of the HI-API-CM continues, software-training efforts will include new solvents, samples, and types of labware, resulting in more robust classification and therefore more faithful artifact assignments.

As previously described and shown in Figure 3, HI-API-CM is able to exploit IR illumination to detect precipitate features, which are otherwise in VIS to the human eye. Combined with the instrument’s high-resolution camera and optics, this is largely responsible for the HI-API-CM’s detection accuracy of 220%; that is, the HI-API-CM routinely identified precipitation artifacts that were missed via manual inspection. In these cases, the precipitation artifacts discovered by the HI-API-CM

were confirmed by manual reinspection of the discrepant wells. In most cases, efforts to manually inspect for precipitate were hampered by a number of challenges including precipitate being occluded by the color of the solvated sample, precipitate, which is a similar color to the labware and precipitate, which had not settled on the bottom of the plate.

HI-API-CM’s fidelity in color classification reveals its sensitivity as a colorimeter. Disagreement in automated versus manual color assignment was not a random error, but rather systematically associated with the presence of faint-colored compound samples. The HI-API-CM performs color analysis in three steps: (1) color detection/resolution, (2) color intensity quantification, and (3) color classification. The ability of these three steps to assign correct values is primarily dependent on the concentration of the colored sample in the well. Wells containing samples with deep color, that is, with intensity above the HI-API-CM’s colorimetric limit of detection (LOD), are detected and classified without any issue. However, as a sample’s concentration is reduced, its concentration approaches the HI-API-CM’s LOD. Once the LOD threshold is exceeded, the HI-API-CM is unable to detect color and therefore the software does not initiate any further color processing. To better understand the HI-API-CM’s




Figure 5. Results of HI-API-CM analysis, as viewed through the UI. An unprocessed image (top left) is presented with other annotated images resulting from HI-API-CM analysis. The UI displays a circle over wells, which HI-API-CM has identified as problematic. An image of wells annotated as containing colored sample (top right), wells containing precipitate image (bottom left), and empty wells image (bottom right) provide a quick reference for users to locate artifacts within the MTP. The UI also allows the user to update the HI-API-CM database.

LOD, the instrument was calibrated using five different DMSO-soluble dyes with colors that span the VIS light spectrum (Table 2).

The values in Table 2 show that the HI-API-CM's LOD varies with color of the dye, with a range of $\sim 10^0$ – $10^{1.6}$ μM . As expected, the fidelity of HI-API-CM color detection above the LOD is 100%. Nevertheless, the 94.0% agreement between the manual and HI-API-CM determinations is acceptable, given that faint-colored samples are less likely to interfere with common HTS-assay detection formats. It is also important to note that the HI-API-CM is able

to distinguish wells with more than one artifact, for example, wells containing colored precipitate (Fig. 3) or partially filled wells containing colored samples. Again, this type of sample does not present a classification problem as long as its color intensity is above the HI-API-CM's corresponding LOD.

CONCLUSION

The HI-API-CM platform represents an important paradigm shift in conventional HTS compound library management practices. A competitive comparison with current QA/QC instrumentation is presented in Table 3. In its current version, HI-API-CM can rapidly identify and classify MTP wells containing colored compound samples, empty, filled, and partially filled wells and the presence of precipitate/crystals in MTP wells with high fidelity. HI-API-CM's rapid, economical, and nondestructive analysis allows the monitoring of an entire HTS compound collection through its lifecycle,

Table 1. Results of manual (visual) and HI-API-CM pilot study ($n = 5760$ wells total)

Artifact type	Number of wells assigned from manual (visual) inspection	Number of wells assigned by HI-API-CM analysis	Agreement, %
Empty	1620	1617	99.8
Filled	4122	4126	100.1 ^a
Partially filled	18	17	94.4
Precipitate	5	11	220.0 ^a
Color	2050	1927	94.0 ^b

^aThe HI-API-CM was found to have greater accuracy than manual methods. Visual reinspection of affected plates confirmed the presence of HI-API-detected artifact.

^bAgreement increases to 100% when the HI-API-CM colorimetric LOD is taken into account.

Table 2. HI-API-CM LOD values for different chromophores

Color	Dye name	LOD (μM)
Red	Allura red	35.1
Green	Brilliant green	9.0
Blue	Euroglucine	0.6
Yellow	Tartrazine	12.7
Violet	Gentian violet	1.3

Table 3. Capabilities of different analytical methods and instrumentation used for quality control of HTS compound libraries

Analytical instrument or method	Precipitate	Empty wells	Presence of bubbles	Sample volume	MW ID	Sample color	Analysis time, 384-well plate
Nephelometry ^{10,11,16}	Yes	No	No	No	No	No	~ 5 min
Acoustic devices ^{12–15}	^a	Yes	^a	Yes	No	No	~ 7 min ^a
LC–MS ⁹	No	No	No	No	Yes	Yes	3–40 h ^b
HI-API-CM	Yes	Yes	Yes	Yes	No	Yes	~ 1 min

^aNot designed to unambiguously detect problem but will flag acoustic anomalies; the effect of suspended solids/colloids are unknown. Acoustic methods require proprietary consumables.

^bAnalysis requires destructive sampling.

MW ID, molecular weight identification.

a practice not currently possible. Future enhancements of the HI-API-CM include the implementation of *chemical imaging* (CI) for quantitative analyses of MTP wells. CI is a promising technology that can readily provide greater quantitative precision of sample volume and DMSO hydration levels in MTPs.

ACKNOWLEDGMENTS

The work of PB, LS, RE, and PH was supported by the Scripps Florida Funding Corporation and by the National Institutes of Health Molecular Library Screening Center Network (MH084512).

Competing Interests Statement: The authors certify that they have no relevant financial interests in this article, and that any/all financial and material support for this research and work are clearly identified in the article.

REFERENCES

- Cox, B.; Denver, J. C.; Binnie, A.; Donnelly, M. C.; Evans, B.; Green, D. V.; Lewis, J. A.; Mander, T. H.; Merrit, A. T.; Valler, M. J.; Watson, S. P. Application of high-throughput screening techniques to drug discovery. *Prog. Med. Chem.* **2000**, *37*, 83–133. PMID: 10845248.
- Yasgar, A.; Shinn, P.; Jadhav, A.; Auld, D.; Michael, S.; Zheng, W.; Austin, C. P.; Inglese, J.; Simeonov, A. Compound management for quantitative high-throughput screening. *JALA Charlottesville, Va.* **2008**, *13*, 79–89.
- Comley, J. Compound management in pursuit of sample integrity. *Drug Discov. World* **2005**, *6*, 59–78.
- Chan, J. A.; Hueso-Rodriguez, J. A. Compound library management. *Methods Mol. Biol.* **2002**, *190*, 117–127.
- Spencer, P. A. The challenges of managing a compound collection. *Laboratory Automation, Eur. Pharm. Rev.* **2004**, *9*, 51–57.
- Kozikowski, B. A.; Burt, T. M.; Tirey, D. A.; Williams, L. E.; Kuzmak, B. R.; Stanton, D. T.; Morand, K. L.; Nelson, S. L. The effect of freeze/thaw cycles on the stability of compounds in DMSO. In *7th Annual Conference of the Society of Biomolecular Screening*, Baltimore, MD, 2001.
- Kozikowski, B. A.; Burt, T. M.; Tirey, D. A.; Williams, L. E.; Kuzmak, B. R.; Stanton, D. T.; Morand, K. L.; Nelson, S. L. The effect of freeze/thaw cycles on the stability of compounds in DMSO. *J. Biomol. Screen.* **2003**, *8*, 210–215.
- Cheng, X.; Hochlowski, J.; Tang, H.; Hepp, D.; Beckner, C.; Kantor, S.; Schmitt, R. Studies on repository compound stability in DMSO under various conditions. *J. Biomol. Screen.* **2003**, *8*, 292–304.
- Cunliffe, J. M.; Adams-Hall, S. B.; Maloney, T. D. Evaluation and comparison of very high pressure liquid chromatography systems for the separation and validation of pharmaceutical compounds. *J. Sep. Sci.* **2007**, *30*(8), 1214–1223. PMID: 17595957.
- Bevan, C. D.; Lloyd, R. S. A high-throughput screening method for the determination of aqueous drug solubility using laser nephelometry in microtiter plates. *Anal. Chem.* **2000**, *72*(8), 1781–1787.
- Goodwin J. J. Flow Cell System for Solubility Testing (Becton, Dickinson and Co., Franklin Lakes, NJ). U.S. Patent 6,788,409 B2, 2004.
- Olechno, J.; Ellson, R.; Browning, B.; Stearns, R.; Mutz, M.; Travis, M.; Qureshi, S.; Shieh, J. Acoustic auditing as a real-time, non-invasive quality control process for both source and assay plates. *Assay. Drug Dev. Technol.* **2005**, *3*, 425–437.
- Huser, J. *High-Throughput Screening in Drug Discovery Methods and Principles in Medicinal Chemistry*. Wiley-VCH InterScience publications: Weinheim, Germany; 2006; pp 37–73.
- Labcyte 1190 Borregas Ave Sunnyvale CA 94089 Echo Qualified Microplate, <http://www.labcyte.com/> (accessed March 23, 2011).
- EDC Biosystems 1804 McCarthy Blvd Milpitas, CA 95035 ATS-100 Acoustic Transfer System, <http://www.edcbiosystems.com> (accessed March 23, 2011).
- IMGEN Technologies 312 S. Washington St, ste 3E Alexandria VA 22314 BMG NEPHELOstar, <http://www.imgen.com/nephelostar.htm> (accessed March 23, 2011).